

## HypaPlast Plasters

### HypaPlast Fabric Plasters




#### Quality fabric plasters in all sizes

- + High quality hypoallergenic plasters
- + Individually wrapped and sterile
- + Comfortable stretch fabric conforms easily to body parts
- + Breathable



### HypaPlast Pink Washproof Plasters




#### Washproof plasters at a great price

- + High quality hypoallergenic plasters
- + Individually wrapped and sterile
- + Washproof



### HypaPlast Blue Catering Plasters



#### Essential for catering environments

- + High quality hypoallergenic plasters
- + Individually wrapped and sterile
- + Washproof and metal detectable
- + Essential for catering and food handling applications




Description	A Fabric	B Washproof	C Blue Metal Detectable	Price €
Finger 7.2x2.5cm Pk of 100	404652	407295	400682	12.20
Knee 7.2x5cm Pk of 100	404655	407299	407962	21.00
Anchor Pk of 100	404102	500034	407903	21.00
Finger Tip Pk of 100	404106	407292	400685	18.50
Assorted Pk of 100	404688	407390	400681	14.95
Finger Extension 2x12cm Pk of 40	-	-	479070	8.90

## HypaPlast Heavy Duty Fabric Plasters

#### Long lasting heavy duty plasters

- + High quality hypoallergenic plasters



Code	Description	Price €
404752	Finger 2.5x7.5cm Pack of 100	16.50
404788	Assorted Pack of 100	17.95

## Coverplast® Plasters



### Coverplast® Classic Fabric Plasters



#### Dependable adhesive wound dressing

- + Good and lasting fixation from a strong rubber-based adhesive
- + Conformable
- + Breathable
- + Low-adherent wound pad
- + Robust and durable
- + Substrate made from 100% cotton




### Coverplast® Barrier Waterproof Plasters



#### Effective infection protection

- + Low-adherent absorbent pad
- + Waterproof
- + Breathable, reducing the risk of skin maceration
- + Wide, all-round adhesive margin
- + Viral and bacterial barrier
- + Breathable polyurethane film
- + Impermeable to oil and fat
- + Skin-friendly acrylic adhesive
- + Highly conformable




### Coverplast® Detectable X-Ray Plasters



#### Easily identifiable and reliable detection

- + X-ray by tungsten filaments
- + Metal detectable (non-ferrous) by aluminium content
- + Substrate made of thin, highly elastic polyurethane film
- + Absorbent, low adherent wound pad
- + Skin-friendly acrylic adhesive
- + Waterproof and breathable
- + Impermeable to oil and fat
- + Viral and bacterial barrier



Description	D Classic Fabric	E Barrier Waterproof	F Blue Metal Detectable	Price €
Finger 2.5x7.5cm Pk of 100	404852	407495	407905	17.50
Knee 5x7.2cm Pk of 100	-	-	407907	24.25
Anchor Pk of 100	404855	-	-	22.50
Assorted Pk of 126	404888	-	-	22.50
Assorted Pk of 120	-	407490	-	22.50
Assorted Pk of 95	-	-	407904	22.50